

465

GÖZ KORUMALI MASKE

Surgeon's face Mask

Mascaras cirurgias

Op-Gesichtsmaske

Mond-en-neusmasker

Mascherina chirurgica

Mascarilla quirurgica

Operations munskydd



Surgeon's face masks with ties

Design: With three plies of non-woven joined together as "sandwich". Each non-woven ply (outer facing, filter media, inner facing) is differently structured in order to satisfy its own specific function; protective, filtering and anti-macerating. These masks are ultrasonically welded.

Classification: Type II according to EN 14683:2005.

Bacteria filtration efficiency= 99,5%

Nose adapter: One long piece ensures a correct face adjustment.

Ties: They are very long.

Color : Green

Packing: Shelf box + Transport carton

Bağcıklı yüz korumalı maske

Dizayn: Değişik fonksiyonları yerine getirmeye yarayan "örgüsüz" yapıda "dört pile" den oluşmuştur. Bu pilelerin yapıları ve işlevleri değişik gösterir ve

- Üst Pile, koruma,
- İkinci Pile, hidropellent-su tutmama
- Üçüncü Pile, filtrasyon,
- Alt Yüz Pilesi ise anti-maserasyon-nem tutmama görevi görür.

Maske ultrasonic yapıştırma yöntemi ile üretilmiş olup, dikişsiz yapıdadır.

Sınıflama: : EN 14683:2005. standartına göre tip II

Bakteri tutma etkinliği = 99,5%

Burun adaptörü: Tek, uzun bir parça maskenin yüze tam olarak ayarlanmasına olanak tanır.

Bağcıklar: 4 uzun bağcığa sahiptir.

Renk : Yeşil

Paketleme :Dispenser Kutu + Taşıma Kutusu

Non Steril / Steril

Kod	Açıklama	Kutu Adeti	Koli Adeti
465-002	Göz koruma kalkanlı	50	600