

4371

DRENAJ-TRAKEOSTOMİ PEDI

Drainage-Tracheotomy sponges

Compressas de drenage-traqueotomia

Compresses de drenage-tracheotomie

Drain-tracheo - kompressen

Splitkompres

Compresa de drenage - traqueotomia

Drainage - trakeal kompress



Composition: Laminate of 3 layers:

-Non adherent net.

-Absorbent fleece.

-Water repellent viscose non-woven.

Approx. weight: 275g/sq.m.

Presentation: Drainage: With an expandable hole to accommodate tubes of diameter from 4,5mm to 9mm.

Tracheotomy: With an expandable hole to accommodate tubes of diameter from 12mm to 16mm.

Packing: Peel-open pack + shelf box + transport carton.

Birleşim: Birbirine "lamine" edilmiş üç teknolojik kattan oluşmaktadır.

o Yaraya yapışmayan mikrogözenekli net,

o Emici özelliğe sahip katman,

o Kan, sekresyon gibi sıvıları dışarı vermeyen,

dokusuz viskon katman.

Yaklaşık ağırlık: 275g / m².

Sunum Şekli: Drenaj pedinde ortasında, çapı 4,5mm – 9mm arasında değişen, **trakeostomi** pedinde ise çapı 12mm – 16mm arasında değişen, değişik ebatlı hortum ve tüpleri yerleştirmeye olanak veren, ayarlanabilir delik bulunmaktadır.

Paketleme: Tekli steril paket + dispansir kutu + taşıma kutusu.

Steril

Kod	Ölçü	Kutu Adeti	Koli Adeti
4371-001	6cm x 7cm	150	2700
4371-002	8cm x 9cm	150	1500

B.V. Bastos Viegas, s.a.

Türkiye yetkili satıcısı: BGM Tıbbi Ürünler-Şemsettin Günaltay Cd. No:73/7 Suadiye İst. - info@bgmegitim.com