

4201

SARGI ALÇI

Plaster of Paris

Bande platree

Gipsbinde

Gipswindsel

Benda gessata

Venda enyesada

Gipsbinda



P.O.P Bandages:

Composition: Leno cotton gauze, without fluorescence, impregnated with gypsum of 88% calcium sulfate

Length: 2 or 3 m.

Approx. weight: 450g/sq.m

Full time to shape: 130s.

Imersion time: 1 to 3s.

Dryaing time: 4 to 7 minutes.

Resistance after 30 minutes: 110Kgf/sq.cm

Packing: Individual wrapping in film + transport carton.

Alçı Özellikleri:

Birleşim: %88 oranında Kalsiyum Sulfat alçıya emdirilmiş, Leno pamuklu bez.

Uzunluk: 2-3 m.

Yaklaşık Ağırlık: 450g / m2

Tam Biçimlendirme Süresi: 130 sn.

Daldırma Süresi: 1-3 sn.

Kuruma Süresi: 4-7 Dakika

30 dakika sonraki direnci: 110kg / cm2

Paketleme: Şeffaf film ambalaj + Taşıma Kutusu

Non Steril

Kod	Ölçü		Koli Adeti
4201-001	2m x 6cm	1	20
4201-002	3m x 10cm	1	20
4201-003	3m x 15cm	1	20
4201-004	3m x 20cm	1	30

BV Bastos Viegas, s.a.

Türkiye yetkili satıcısı: BGM Medikal-Şemsettin Günaltay Cad. No:61/7 Suadiye İst. - info@bgmegitim.com