

# PHIL-SPLINT™

Option 1. 8 layers of knitted Fiberglass/Polyester in core  
 Option 2. 1 layer of non-woven fabric for in core Coated by polyurethane.



CODE NO. **PS-8XXX**

**PS-1XXX**

Initial for 'PHIL-SPLINT'

Numbers of in-core materials

Description for size of spec.

ORDER NO.	APPLICATION	SIZE	SPEC.	PACKING
PS-8210	Short Arm	2'x10'	5.0cmx25.0cm	7EA/INBOX
PS-8312	Short Arm	3'x12'	7.5cmx30.0cm	5EA/INBOX
PS-8415	Short Arm	4'x15'	7.5cmx87.5cm	5EA/INBOX
PS8430	Long Arm/Short Arm	4'x30'	10.0cmx37.5cm	5EA/INBOX
PS8530	Short Arm	5'x30'	10.0cmx75.0cm	5EA/INBOX
PS-8545	Long Leg	5'x45'	12.5cmx75.0cm	5EA/INBOX
PS-8645	Long Leg	6'x45'	12.5cmx112.5cm	5EA/INBOX



# PHIL ROLL SPLINT™

(ORTHOPEdic CASTING TAPE)



## CHARACTERISTICS

- Much more rigid and light
- Wide ranges of light or severe damaged bone-fracture
- Easy application and patients feel comfortable
- Various size (8 size of ready made)




- Various size (5 size)
- No loss to you! Free size on your decision!



International certification marked

Roll clips by sizes



## VARIOUS COLOR OF PACKAGE

Pouch Color - 2": Purple / 3": Olive / 4": Sky Blue / 5": Dark Blue / 6": Deep Red : Color difference of pouch can make you to distinguish the stock of each size at a single glance and pouch & in box color is same for consistent design acknowledgement.

ORDER NO.	SIZE	SPEC.	PACKING
PR-802	2'x180'	5.0cmx450.0cm	1EA/INBOX
PR-803	3'x180'	7.5cmx450.0cm	1EA/INBOX
PR-804	4'x180'	10.0cmx450.0cm	1EA/INBOX
PR-805	5'x180'	12.5cmx450.0cm	1EA/INBOX
PR-806	6'x180'	15.0cmx450.0cm	1EA/INBOX